

For Research Use Only.  
Not For Use In Diagnostics.

# CoraLite488-conjugated Recombinant Rabbit Anti-Mouse IgG Kappa Light Chain



Catalog Number: SA00014-3

## Information

Catalog Number:

SA00014-3

Host:

Rabbit

Applications:

IF

Reactivity:

Mouse

Physical State:

Liquid

Conjugation:

CoraLite®488

## Recommended Dilutions

1:100-1:500

## Preservative

Proclin300

## Fluorophore

CL488 Fluorescent Dye, Amax=488 nm; Emax=515 nm

## Stabilizer

0.5% BSA

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

## Purity

Protein A purification

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

