

For Research Use Only.
Not For Use In Diagnostics.

CoraLite488-conjugated Mouse Anti-Heavy Chain of Rabbit IgG

Catalog Number:SA00014-4



Information

Catalog Number:

SA00014-4

Host:

Mouse

Applications:

IF, FC

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

CoraLite®488

Recommended Dilutions

1:100-1:500

Preservative

Proclin300

Fluorophore

CL488 Fluorescent Dye, Amax=488 nm; Emax=515 nm

Stabilizer

0.5% BSA

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

Purity

Protein G purification

For technical support and original validation data for this product please
contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

