

For Research Use Only.
Not For Use In Diagnostics.

CoraLite Plus 750-conjugated Goat anti-rabbit IgG(H+L)

Catalog Number:SA00015-1



www.ptglab.com

Information

Catalog Number:

SA00015-1

Host:

Goat

Applications:

WB

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

CoraLite® Plus 750

Recommended Dilutions

1:5000 – 1:10000

Preservative

0.1% Proclin-300

Fluorophore

CoraLite Plus 750, Amax=755±4 nm; Emax=780±4 nm

Stabilizer

10 mg/mL BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 10 mg/mL BSA, 0.1% Proclin-300, pH 7.4.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from culture media supernatant by immunoaffinity chromatography using Protein G beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data