

For Research Use Only

# Rabbit anti-Goat IgG serum

Catalog Number: B900710



[www.ptgcn.com](http://www.ptgcn.com)

---

Physical State	Liquid
Source	Collected from Japanese White Rabbit which has been immunized by Goat IgG. Aseptically collected and processed to minimize hemolytic activity.
Size	1ml
Preservative	0.02% Sodium Azide
Storage	-20°C, avoid repeated thawing and freezing
Note	This product is for research use only, not for diagnostic or therapeutic use.

---

For technical support and original validation data for this product please contact

T: 027-87531629

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.