

For Research Use Only

anti-Mouse IgG3 VHH Agarose for Immunoprecipitation



www.ptgcn.com

Catalog Number: mIG3a

Basic Information

Catalog Number:
mIG3a

Applications:
IP, Co-IP

Conjugate:
Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)

Host:
Alpaca

Type:
Nanobody

Class:
Recombinant - Animal free production

Description

anti-Mouse IgG3 IP Beads is an affinity resin for IP of Mouse IgG3. It consists of rabbit specific VHHs (Nanobodies) coupled to Agarose beads.

Binding capacity

/

Elution buffer

SDS Sample Buffer

Wash buffer compatibility

/

Affinity (K_D)

/

Storage

Storage:
+4°C / do not freeze!

Storage Buffer:
20% Ethanol

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com
W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.