

For Research Use Only

# anti-Mouse IgG Magnetic VHH Agarose for Immunoprecipitation



[www.ptgcn.com](http://www.ptgcn.com)

Catalog Number: mIGma

## Basic Information

**Catalog Number:**

mIGma

**Applications:**

IP, Co-IP

**Conjugate:**

Magnetic Agarose beads; ~40 um (cross-linked 6% magnetic agarose beads)

**Host:**

Alpaca

**Type:**

Nanobody

**Class:**

Recombinant - Animal free production

## Description

anti-Mouse IgG IP Beads is an affinity resin for IP of all subtypes of Mouse IgG. It consists of rabbit specific VHHs (Nanobodies) coupled to Magnetic Agarose beads.

## Binding capacity

/

## Elution buffer

SDS Sample Buffer

## Wash buffer compatibility

/

## Affinity ( $K_D$ )

/

## Storage

**Storage:**

+4°C / do not freeze!

**Storage Buffer:**

20% Ethanol

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [www.ptgcn.com](http://www.ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.