

For Research Use Only

anti-Halo recombinant VHH, unconjugated



www.ptgcn.com

Catalog Number: ot

1 Publications

Basic Information

Catalog Number:

ot

Applications:

Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

Host:

Alpaca

Conjugate:

Unconjugated

Type :

Nanobody

Class:

Recombinant

RRID:

AB_2827559

Molecular Weight:

14.7 kDa

Description

Alpaca anti-Halo VHH, purified recombinant binding protein

Specificity/Target

Halo-tag

Affinity (K_D)

Dissociation constant K_D of 2 nM

Storage

Storage:

Shipped at ambient temperature. Upon receipt store at 4°C; stable for one year. Do not freeze.

Storage Buffer:

25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Structure of Halo.