

For Research Use Only

SARS-CoV-2 Spike Recombinant VHH [NM1266]



www.ptgcn.com

Catalog Number: sc-NM1266

Basic Information

Catalog Number:

sc-NM1266

Clone No.:

NM1266 (NM1230-NM1228 fusion) [Wagner et al. EMBO Rep. 2021 May 5;22(5):e52325. PMID: 33904225]

Applications:

BLI, Multiplex ACE2 competition assay

Host:

Alpaca, recombinantly produced

Conjugate:

Unconjugated

Type:

Nanobody / VHH, biparatopic

Class:

Recombinant

RRID:

AB_2892256

Purification Method:

Recombinant expression, affinity purification by IMAC.
Note: SARS-CoV-2 Spike Recombinant VHH carries a C-terminal His-tag.

Description

Recombinant bivalent Nanobody against the Receptor-binding domain (RBD) of SARS-CoV-2 Spike protein

Affinity (K_D)

NM1230: 9.5 nM
NM1228: 0.75 nM

Storage

Storage:

Aliquot upon receipt and store at -20°C/-4°F. Avoid freeze-thaw cycles.

Storage Buffer:

PBS, 0.09% sodium azide
Safety datasheet (SDS): sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



BLI binding kinetics of SARS-CoV-2 Spike Recombinant VHH [NM1266] to RBD. Biotinylated RBD was immobilized on FortéBio Streptavidin (SA) Biosensors and assayed with 25 nM of SARS-CoV-2 Spike Recombinant VHH [NM1266] (ChromoTek sc-NM1266). For comparison the binding kinetics of 25 nM SARS-CoV-2 Spike Recombinant VHH [NM1228] (ChromoTek sc-NM1228) and SARS-CoV-2 Spike Recombinant VHH [NM1230] (ChromoTek sc-NM1230) are shown.