

For Research Use Only

anti-GST recombinant VHH, unconjugated



www.ptgcn.com

Catalog Number: st

4 Publications

Basic Information

Catalog Number:
st

Applications:
Suitable for conjugation of dyes, biotin, beads,
surfaces etc. via NHS ester reaction

Host:
Alpaca

Conjugate:
Unconjugated

Type :
Nanobody

Class:
Recombinant

RRID:
AB_2631368

Molecular Weight:
14.7 kDa

Description

Alpaca anti-GST VHH, purified recombinant binding protein

Specificity/Target

Glutathione S-Transferase (GST) from *Schistosoma japonicum*

Affinity (K_D)

Dissociation constant K_D of 1 nM

Storage

Storage:

Shipped at room temperature. Upon receipt store at 4°C; stable for one year. Do not freeze.

Storage Buffer:

25mM TAPS pH 8.5, 500 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide


For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com
W: www.ptgcn.com

This product is exclusively available under
Proteintech Group brand and is not available
to purchase from any other manufacturer.

Selected Validation Data



Structure of GST.