

For Research Use Only

anti-V5 recombinant VHH, unconjugated



www.ptgcn.com

Catalog Number: v5t

Basic Information

Catalog Number:
v5t

Applications:
Suitable for conjugation of dyes, biotin, beads, surfaces etc. via NHS ester reaction

Host:
Alpaca

Conjugate:
Unconjugated

Type :
Nanobody

Class:
Recombinant

RRID:
AB_2868499

Molecular Weight:
14.6 kDa

Description

Alpaca anti-V5 VHH, purified recombinant binding protein

Specificity/Target

V5-tag sequence GKIPNPLLGLDST

Affinity (K_D)

Dissociation constant K_D of 40 nM

Storage

Storage:

Shipped at ambient temperature . Upon receipt store at 4°C; stable for one year.

Storage Buffer:

25mM TAPS pH 8.5, 137 mM NaCl, 5 mM EDTA, preservative: 0.09 % sodium azide

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com
W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Structure of VHH.