

For Research Use Only

Vimentin RecAb, VHH-rabbit IgG Fc fusion



www.ptgcn.com

Catalog Number: **vfrb**

Basic Information

Catalog Number:

vfrb

Clone No.:

CTK0211

Applications:

IF

Host/IsoType:

Rabbit/IgG

Conjugate:

Unconjugated

Type:

Primary antibody; chimeric fusion of Nanobody (VHH) and rabbit IgG Fc domain

Class:

Recombinant

RRID:

AB_2889373

Purification Method:

Recombinantly expressed and purified by IMAC

Description

Rabbit monoclonal antibody to Vimentin. Recombinant fusion of anti-Vimentin VHH [CTK0211] and rabbit IgG Fc domain.

Specificity/Target

Vimentin

Recommended Dilution

IF/ICC: 1:400 - 1:2,000

Storage

Storage:

Shipped on ice. Aliquot upon receipt and store at -20°C/-4°F. Avoid freeze-thaw cycles.

Storage Buffer:

PBS, 10% glycerol, Preservative: 0.02 % sodium azide, safety datasheet (SDS): sodium azide

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MDCK cells were immunostained with Vimentin recombinant antibody, VHH-rabbit IgG Fc fusion [CTK0211] (vfrb, 1:1,000) and Nano-Secondary? alpaca anti-human IgG/anti-rabbit IgG, recombinant VHH, Alexa Fluor? 647 [CTK0101, CTK0102] (srbAF647-1, 1:1,000). Scale bar, 10 μ M.